
CEF4 Peptide

Subcategory: Synthetic Peptide

Cat. No.: 350123

Unit: 2 mg

Description:

CEF4 corresponds to aa 342-351 of the influenza A virus nucleocapsid protein.

Format: Each vial contains 1 mg of lyophilized solid packaged under an inert gas and supplied as a trifluoroacetate salt.

Alternate Names: Nucleocapsid protein; nucleoprotein; protein N; NP

Accession No.: P31609

MW: 1148.42 g/mol

Sequence: H-Arg-Val-Leu-Ser-Phe-Ile-Lys-Gly-Thr-Lys-OH or H-RVLSFIKGTK-OH

Composition: C53H93N15O13

Purity: Purity > 95% by HPLC

Solubility: Distilled water for a solution up to 2 mg/ml, otherwise we recommend using acetonitrile.

Storage: Store at -20°C. The product is hygroscopic and must be protected from light. Product is guaranteed one year from the date of shipment. Following reconstitution, aliquot and store at -20°C.

For research use only, not for diagnostic or therapeutic procedures.